



www.protea-vital.com


The expectations of the final beneficiaries are the basis of our entire activity. By understanding them very well, we strive through honest work and dedication to fulfill and to overcome their challenges constantly. In our work we respect all applicable quality standards, and the requirements of our customers are sacred to us. We are constantly working to offer the best solutions, and the processes through which we achieve this is ultimately leading us to be a robust, reliable, world class business.




PROTEA VITAL

Designed to provide a person with a stoma with improved comfort, quality of life, increased privacy and a lifestyle as close as possible to that before surgery.

Our products are carefully made, to your high expectations and quality standards.
ISO 13485 certified





PROTEA VITAL CF3L

The CF3L is a unitary (one piece), closed end (no evacuation, closed) ostomy collection bag with carbon filter. All components of the product are removed from the body once the bag is changed. This ostomy collection bag does not allow effluents to be discharged, followed by prolonged use.

Learn more about PROTEA VITAL CF3L's product specifications at www.protea-vital.com/CF3L



PROTEA VITAL DCF3L

The PROTEA VITAL DCF3L is a unitary (one piece) ostomy collection bag with emptying (with discharge, open-end or drainable) – using a clamp closure. A high-capacity bar filter is used for gas expenditure. All components of the product are removed from the body once the bag is changed. This ostomy collection bag allows the discharge of effluents by means of a special portion for emptying, which closes / opens with a clamp.

Learn more about PROTEA VITAL DCF3L's product specifications at www.protea-vital.com/DCF3L



PROTEA VITAL DVF3L

The PROTEA VITAL DVF3L is a unitary (one piece) ostomy collection bag with emptying (with discharge, open-end or drainable) – using a Velcro closure. A high-capacity bar filter is used for gas expenditure. All components of the product are removed from the body once the bag is changed. This ostomy collection bag allows the discharge of effluents by means of a special portion for emptying, which closes/opens with a Velcro system.

Learn more about PROTEA VITAL DVF3L's product specifications at www.protea-vital.com/DVF3L